



ALS Bottle Kits Save Time and Money

Prepared with the client and sample collection personnel in mind



ALS prepares and delivers bottle kits with one goal in mind, to make the sample collection process as efficient and problem-free as possible.

Unless specified otherwise, sample containers are pre-labeled and pre-preserved. Sample bottles that contain a chemical preservative include a sticker on the container that identifies the preservative. The sample labels that accompany the sample kits include spaces for client identification, sample identification, date and time of collection, sample collector identification, analysis(es) to be performed, and preservative identification. Glass containers are enclosed in bubble bags to ensure container integrity. VOC trip blanks and cooler temperature blanks are clearly labeled for easy identification.

Sample containers are typically grouped and packaged such that all of the containers required for a given sample collection point are enclosed in a single zippered storage bag (or multiple bags if the number of sample containers requires it). This eliminates the need to pull individual containers from multiple coolers at each sample collection point. Cooler liners are provided for packing and icing of samples for transport to the laboratory, and chain of custody documents are included for field use. Custody seals are also provided for cooler and/or sample container use.

Our bottle kits include:



Chemical preservative stickers and sample labels with each parameter clearly identified.



Pre-labeled containers are grouped by sample location and enclosed in a zippered storage bag.



Clean, lined coolers neatly packaged for delivery to clients. Sample locations are clearly marked for quick access when sampling.



Pre-printed chain of custody forms with client and project names to minimize time spent in the field.



After a final quality control check, each cooler is sealed with our Quality Seal of Excellence identifying the client and project prior to shipping.

ALS delivers the best, and we are proud of our stellar track record. We will not settle for the status quo when performance excellence is our goal. ALS recognizes you have a choice, and we are committed to your success.

Call today to schedule a laboratory tour and speak with one of our Project Managers to learn how we can assist you on your testing requirements.



10450 Stancliff Road, Suite 210
Houston, TX 77099
+1 281 530 5656